



## **Troubleshooting Guide – *COSMO* cDNA synthesis Kit *COSMO* RT-PCR Master Mix**

The following guide can be used to troubleshoot any PCR reactions using the *COSMO* cDNA mix or *COSMO* RT-PCR Master Mix.

If you require any additional information, please email us at [info@willowfort.co.uk](mailto:info@willowfort.co.uk).

<u>Problem</u>	<u>Causes</u>	<u>Possible Actions</u>
<b>No or low Amplification</b>	Poor quality RNA	<ul style="list-style-type: none"> <li>• Reduce the freeze-thaw cycles, to reduce degradation</li> <li>• Check the RIN score before cDNA synthesis</li> <li>• Ensure a clean working area to minimise RNase contamination</li> <li>• Use nuclease-free water</li> </ul>
	Poor RNA purity	<ul style="list-style-type: none"> <li>• Consider sample specific purification kits for RNA</li> <li>• Measure RNA purity using a UV spectrophotometer.</li> <li>• Reduce the concentration of inhibitors by dilute RNA in nuclease free water</li> <li>• Remove salts and inhibitors by repeating purification.</li> </ul>
	Low quantity of RNA	<ul style="list-style-type: none"> <li>• Check RNA concentrations using a UV spectrophotometer.</li> <li>• Consider using reagents and reverse-transcriptase with increased sensitivity to low RNA concentrations.</li> </ul>
	High GC content	<ul style="list-style-type: none"> <li>• Use thermostable reverse transcriptase to allow for higher temperatures, breaking the GC bonds.</li> </ul>
	High secondary structure content	<ul style="list-style-type: none"> <li>• Increase temperature during reverse transcription to prevent hairpin structures.</li> <li>• Before reverse transcription, prevent secondary structures by incubating RNA at a high temperature, followed by a transfer to ice.</li> </ul>
	Insufficient reverse transcriptase	<ul style="list-style-type: none"> <li>• Increase cDNA quantity</li> <li>• Consider a reverse transcriptase with increased sensitivity</li> <li>• Consider a reverse transcriptase efficient against potential inhibitors</li> </ul>

<u>Problem</u>	<u>Causes</u>	<u>Possible Actions</u>
----------------	---------------	-------------------------

**No or low Amplification continued...**

Insufficient time and temperature for optimum Reverse Transcriptase reaction

- Check RT-PCR machine for correct time and temperature

Error in primer design

- Choose primers specific to RNA i.e. if using degraded RNA, random primers would be beneficial rather than oligo(dT) primers

Reagent quality

- Refer to manufacturer's storage suggestions
- Only use fresh reagents with reverse transcriptase
- Reduce the number of freeze-thaw cycles
- Ensure a complete mix of reagents
- Check for any precipitation in any reagents.

DNA contamination

- Use a control e.g. No-RT control
- Treat RNA samples before beginning reverse transcription

**Non-Specific Amplification**

Error with primer design

- Increase primer length to increase specific amplification.
- Increase annealing temperature to improve specificity of primer binding.

Poor RNA quality

- Check RNA concentration before cDNA synthesis
- Reduce the number of freeze-thaw cycles
- Include a RNase inhibitor during reverse transcription
- Use nuclease-free water buffer when storing RNA
- Include a genomic DNA step during RNA extraction

**Truncated cDNA**

High secondary structure content

- Increase temperature during reverse transcription to prevent hairpin structures.
- Before reverse transcription, prevent secondary structures by incubating RNA at a high temperature, followed by a transfer to ice.

Errors with primers

- Choose primers specific to RNA i.e. if using degraded RNA, random primers would be beneficial rather than oligo(dT) primers
- Use oligo(dT) when synthesizing full-length cDNA
- Optimise primer concentration

<u>Problem</u>	<u>Causes</u>	<u>Possible Actions</u>
<b><i>PCR Reaction Not Reproducible or Reaction Stopped Working</i></b>	Different cycling conditions	<ul style="list-style-type: none"> <li>• Use the same thermal cycler for optimization and all future experiments.</li> </ul>
	dNTPs, primers or sample degraded	<ul style="list-style-type: none"> <li>• Aliquot components to reduce number of freeze-thaw cycle.</li> </ul>
	Error in set up	<ul style="list-style-type: none"> <li>• Ensure correct reagents have been added</li> <li>• Check correct thermal cycler program has been selected.</li> </ul>
	Change in component	<ul style="list-style-type: none"> <li>• Check for new components added to the master mix (e.g. new batch of primers)</li> </ul>
	Inhibitors in template	<ul style="list-style-type: none"> <li>• Decrease template concentration.</li> </ul>