



Troubleshooting Guide –

HERA RT-qPCR Kit

HERA SYBR[®] Green RT-qPCR kit

The following guide can be used to troubleshoot any reactions using the *HERA RT-qPCR Kit* or *HERA SYBR[®] Green RT-qPCR kit*.

If you require any additional information, please email us at info@willowfort.co.uk.

<u>Problem</u>	<u>Causes</u>	<u>Possible Actions</u>
No or low Amplification	Poor quality RNA	<ul style="list-style-type: none"> • Reduce the freeze-thaw cycles, to reduce degradation • Check the RIN score before cDNA synthesis • Ensure a clean working area to minimise RNase contamination • Use nuclease-free water
	Poor RNA purity	<ul style="list-style-type: none"> • Consider sample specific purification kits for RNA • Measure RNA purity using a UV spectrophotometer. • Reduce the concentration of inhibitors by dilute RNA in nuclease free water • Remove salts and inhibitors by repeating purification.
	Low quantity of RNA	<ul style="list-style-type: none"> • Check RNA concentrations using a UV spectrophotometer. • Consider using reagents and reverse-transcriptase with increased sensitivity to low RNA concentrations.
	High GC content	<ul style="list-style-type: none"> • Use thermostable reverse transcriptase to allow for higher temperatures, breaking the GC bonds.
	High secondary structure content	<ul style="list-style-type: none"> • Increase temperature during reverse transcription to prevent hairpin structures. • Before reverse transcription, prevent secondary structures by incubating RNA at a high temperature, followed by a transfer to ice.
	Insufficient reverse transcriptase	<ul style="list-style-type: none"> • Increase cDNA quantity • Consider a reverse transcriptase with increased sensitivity • Consider a reverse transcriptase efficient against potential inhibitors

<u>Problem</u>	<u>Causes</u>	<u>Possible Actions</u>
No or low Amplification continued...	Insufficient time and temperature for optimum Reverse Transcriptase reaction	<ul style="list-style-type: none"> • Check RT-PCR machine for correct time and temperature
	Error in primer design	<ul style="list-style-type: none"> • Choose primers specific to RNA i.e. if using degraded RNA, random primers would be beneficial rather than oligo(dT) primers
	Reagent quality	<ul style="list-style-type: none"> • Refer to manufacturer's storage suggestions • Only use fresh reagents with reverse transcriptase • Reduce the number of freeze-thaw cycles • Ensure a complete mix of reagents • Check for any precipitation in any reagents.
Non-Specific Amplification	DNA contamination	<ul style="list-style-type: none"> • Use a control e.g. No-RT control • Treat RNA samples before beginning reverse transcription
	Error with primer design	<ul style="list-style-type: none"> • Increase primer length to increase specific amplification. • Increase annealing temperature to improve specificity of primer binding.
Truncated cDNA	Poor RNA quality	<ul style="list-style-type: none"> • Check RNA concentration before cDNA synthesis • Reduce the number of freeze-thaw cycles • Include a RNase inhibitor during reverse transcription • Use nuclease-free water buffer when storing RNA • Include a genomic DNA step during RNA extraction
	High secondary structure content	<ul style="list-style-type: none"> • Increase temperature during reverse transcription to prevent hairpin structures. • Before reverse transcription, prevent secondary structures by incubating RNA at a high temperature, followed by a transfer to ice.
	Errors with primers	<ul style="list-style-type: none"> • Choose primers specific to RNA i.e. if using degraded RNA, random primers would be beneficial rather than oligo(dT) primers • Use oligo(dT) when synthesizing full-length cDNA • Optimise primer concentration

<u>Problem</u>	<u>Causes</u>	<u>Possible Actions</u>
----------------	---------------	-------------------------

Higher Ct value than expected

RNA template is damaged or degraded

- Review RNA extraction protocol

RNase contamination

- Ensure a clean work area. Use 70% ethanol and dH₂O to clean all equipment and space.
- Add a RNase inhibitor after cDNA synthesis

RT inhibition present

- Use an additional wash of ethanol to remove inhibitors during RNA extraction.

Poor cDNA synthesis

- Check correct temperature on Thermocycler
- If necessary, increase cDNA synthesis temperatures

Low Ct value than expected or lower Ct value than positive control

Template contamination

- Identify source of contamination by checking and replacing reagents.
- Use separate pipettes for pre- and post-cDNA reactions
- Ensure a clean work area. Use 70% ethanol and dH₂O to clean all equipment and space

Concentration of RNA sample are too high

- Confirm the RNA concentration.
- Decrease concentrations of RNA in PCR set-up

Different cycling conditions

- Use the same thermal cycler for optimization and all future experiments.

PCR Reaction Not Reproducible or Reaction Stopped Working

dNTPs, primers or sample degraded

- Aliquot components to reduce number of freeze-thaw cycle.

Error in set up

- Ensure correct reagents have been added
- Check correct thermal cycler program has been selected.

Change in component

- Check for new components added to the master mix (e.g. new batch of primers)

Inhibitors in template

- Decrease template concentration.