



## **Troubleshooting Guide –**

***NOVA* Blood DNA Extraction kit**

***NOVA* Tissue and Mammalian Cell**

**Culture DNA Extraction Kit**

***NOVA* Viral RNA and DNA Extraction kit**

The following guide can be used to troubleshoot any DNA extractions using the *NOVA* kits.

If you require any additional information, please email us at [info@willowfort.co.uk](mailto:info@willowfort.co.uk).

<u>Problem</u>	<u>Causes</u>	<u>Possible Actions</u>
<b><i>DNA is poor quality</i></b>	Starting material was not stored correctly	<ul style="list-style-type: none"> <li>• Nucleases may have degraded your DNA sample.</li> <li>• Store all DNA samples at -20°C or less.</li> <li>• Do not store DNA in water for long periods of time. Instead use TE buffer.</li> <li>• DNA may have been fragmented, making it impossible to use for studies requiring high-molecular-weight DNA applications.</li> </ul>
	Genomic DNA was not stored correctly	
	DNA was extracted from fixed tissue	
<b><i>DNA yield is poor</i></b>	Starting sample was too small	<ul style="list-style-type: none"> <li>• Try using a large volume/weight of sample.</li> <li>• If a large volume/weight is not possible, amplification may be required.</li> <li>• Nucleases may have degraded your DNA sample.</li> </ul>
	Cells were not lysed thoroughly	
	Poor DNA extraction	
<b><i>DNA sample is cloudy</i></b>	Transfer of cellular debris with DNA sample	<ul style="list-style-type: none"> <li>• Centrifuge sample again, taking care to remove only the DNA sample</li> </ul>
	Cellular debris not separate from DNA sample	
<b><i>The A<sub>260</sub>/A<sub>280</sub> ratio does fall within the recommended range</i></b>	A <sub>260</sub> /A <sub>280</sub> is < 1.7, therefore proteins or chaotropic agents are present.	<ul style="list-style-type: none"> <li>• Avoid overloading the silica-based columns.</li> <li>• If your sample has high concentrations of protein, treatment of proteases is recommended before purification.</li> <li>• Follow washing steps closely to prevent carryover of contaminants</li> </ul>
	A <sub>260</sub> /A <sub>280</sub> is > 1.9, therefore RNA could be present.	

<u>Problem</u>	<u>Causes</u>	<u>Possible Actions</u>
<b><i>Problems with downstream applications</i></b>	Poor quality DNA	<ul style="list-style-type: none"> <li>• Nucleases may have degraded your DNA sample.</li> <li>• Store all DNA samples at -20°C or less.</li> <li>• Do not store DNA in water for long periods of time. Instead use TE buffer.</li> <li>• DNA may have been fragmented, making it impossible to use for studies requiring high-molecular-weight DNA applications.</li> </ul>
	Presence of contaminants	<ul style="list-style-type: none"> <li>• Avoid overloading the silica-based columns.</li> <li>• If your sample has high concentrations of protein, treatment of proteases is recommended before purification.</li> <li>• Follow washing steps closely to prevent carryover of contaminants</li> <li>• Repeat DNA precipitation with salts and ethanol.</li> </ul>
	Presence of residual ethanol	<ul style="list-style-type: none"> <li>• Ensure silica filters have been thoroughly dried before eluting sample.</li> </ul>