



Troubleshooting Guide –

NOVA Blood RNA Extraction kit

***NOVA Tissue and Mammalian Cell Culture
RNA Extraction Kit***

NOVA Wilozole RNA Extraction kit

**The following guide can be used to troubleshoot any
RNA extractions using the *NOVA* kits.**

If you require any additional information, please email us at
info@willowfort.co.uk.

<u>Problem</u>	<u>Causes</u>	<u>Possible Actions</u>
Genomic DNA within RNA sample	Insufficient degradation of DNA during lysis	<ul style="list-style-type: none"> • Treat samples with DNase, but make sure DNase is removed after the reaction is complete
	Silica spin columns overloaded	<ul style="list-style-type: none"> • Follow the sample size suggestions to prevent overload
	Starting sample was too small	<ul style="list-style-type: none"> • Use the suggested starting material • Do not overload the spin columns
	Cells were not lysed thoroughly	<ul style="list-style-type: none"> • Optimise the lysis step, particularly samples with tough structures
Low RNA yield	Sample was not stored correctly	<ul style="list-style-type: none"> • Endogenous RNase may have degraded RNA. Adding chaotropic agents will inactivate any RNase present.
	RNase degraded RNA sample during or after extraction	<ul style="list-style-type: none"> • Ensure the correct protective clothing are used during extraction • Make sure consumables are RNase free • Prepare an RNase free environment by cleaning surround area with 70% ethanol and dH₂O
	Sample was not stored correctly	<ul style="list-style-type: none"> • Endogenous RNase may have degraded RNA. Adding chaotropic agents will inactivate any RNase present.
Low RNA integrity	RNase degraded RNA sample during or after extraction	<ul style="list-style-type: none"> • Ensure the correct protective clothing are used during extraction • Make sure consumables are RNase free • Prepare an RNase free environment by cleaning surround area with 70% ethanol and dH₂O
	Extracted total RNA was not stored correctly	<ul style="list-style-type: none"> • It is advised that due to its instability, RNA should be stored at -80°C to reduce risk of degradation.

<u>Problem</u>	<u>Causes</u>	<u>Possible Actions</u>
----------------	---------------	-------------------------

Inhibitors in RNA sample

	A_{260}/A_{280} is < 1.7 , therefore proteins or chaotropic agents are present.	<ul style="list-style-type: none"> • Avoid overloading the silica-based columns. • If your sample has high concentrations of protein, treatment of proteases is recommended before purification. • Follow washing steps closely to prevent carryover of contaminants
	A_{260}/A_{230} is < 1.0	<ul style="list-style-type: none"> • Follow washing steps closely to reduce risk of salt carryover. Increase number of washes if necessary. • Wash sample in ethanol to remove salts.
	Poor quality RNA	<ul style="list-style-type: none"> • RNases may have degraded your DNA sample. • Store all RNA samples at -80°C or less.

Problems with downstream applications

	Presence of contaminants	<ul style="list-style-type: none"> • Avoid overloading the silica-based columns. • If your sample has high concentrations of protein, treatment of proteases is recommended before purification. • Follow washing steps closely to prevent carryover of contaminants • Repeat RNA precipitation with salts and ethanol.
	Presence of residual ethanol	<ul style="list-style-type: none"> • Ensure silica filters have been thoroughly dried before eluting sample.